

Lite



Rediscover your freedom...

Ultra-lightweight, folding
travel scooter



Scan for more
information



eFOLDi



...and travel with confidence.

The eFOLDi Lite is the ultimate solution for inclusive independence

As a recognised leader in the design and manufacture of portable scooters and powered wheelchairs, eFOLDi strives to help people rediscover their freedom and independence.

eFOLDi's award-winning mobility solutions are not only easy to store and transport, they help make travel simple and achievable – whether that's by air, land or sea.



eFOLDi



Fold and go further by car, bus, train, plane, or cruise liner - it's your choice

The eFOLDi Lite is one of the world's lightest, most compact folding travel-scooters.

Weighing just 14kg*, the eFOLDi Lite folds and unfolds in just a few seconds so is easy to store and transport. No complex assembly is required - just unfold and go.

What makes eFOLDi unique...

Lightweight – carry and lift with ease



Compact – folds down to cabin luggage size



Portable – easy to take with you, anywhere



Robust – durable materials give peace of mind



Lite



Compact, comfortable, quick and simple to fold... the Lite combines it all

- Easy to use, folding travel scooter
- Exceptionally lightweight at only 14kg*, yet remarkably robust carrying passengers up to 120kg/19 Stone
- Nimble and very manoeuvrable, with auto turning speed controller for extra stability
- Folds with ease, fits into most car boots and is Airline-safe
- The perfect travel companion!



The lighter way to empower your life



Weighs only 14kg*, for easy storage and transportation

"My eFOLDi Lite folds up like a cabin-sized suitcase in just a few seconds so it can be neatly stored away either in the house or when traveling by car, train, bus or coach."

Roger, Essex



Lite with red seat cushions and arm-rests (optional)

*Excluding battery

Shown with arm-rests (optional)

Simple and intuitive controls

Highly manoeuvrable

Robust, carries users up to 120kg/265lbs

Ultra-lightweight at only 14kg/31lbs*



4-point floating seat suspension system

Battery charging point

Auto turning speed controller

Lithium-ion battery with 10 mile/16km range

Quick release wheels saves another 1kg



Patented, award-winning design

Folds for easy storage or transportation in seconds



Specifications*

Maximum speed	4mph / 6kph
Maximum range	10 miles / 16km
Maximum user weight	19 Stone / 120kg / 265lbs
Climbing ability	1:10 gradient
Folded dimensions	26 x 22 x 15 inches / 65 x 55 x 37cm
Weight (excl. battery)	31lbs / 14kg
Weight (incl. battery)	35lbs / 16kg
Type of use	Pavement only
Braking system	Electro-magnetic
Turning speed controller	50% speed reduction at 20 degree turn or greater
Airline friendly	Yes

* Specifications may vary due to manufacturing tolerances

Lite

...and enjoy being you.

"What a fantastic device. It's been so amazing for my husband to have something so portable"

Lisa Ruxton



Scan for more information

0800 999 8123

sales@efoldi.com www.efoldi.com

Your local eFOLDi retailer:



SunTech UK Limited t/a eFOLDi, 25 Ormside Way, Holmethorpe Industrial Estate, Redhill, Surrey RH1 2LW.

eFOLDi