

Power Chair



Rediscover your freedom...

The ultra-lightweight folding powered wheelchair



Scan to see the Power Chair in action



eFOLDi



...and travel with confidence.

The eFOLDi Power Chair is the ultimate solution for inclusive independence

As a recognised leader in the design and manufacture of portable scooters and powered wheelchairs, eFOLDi strives to help people rediscover their freedom and independence. eFOLDi's award-winning mobility solutions are not only easy to store and transport, they help make travel simple and achievable – whether that's by air, land or sea.



eFOLDi



Fold and go further by car, bus, train, plane or cruise liner - it's your choice

The eFOLDi Power Chair is one of the world's lightest, most compact folding powered wheelchairs. Weighing just 15kg*, the ultra-portable Power Chair folds and unfolds in just a few seconds so is easy to store and transport, making the perfect travel companion. No complex assembly is required - just unfold and go.

What makes eFOLDi unique...

Lightweight – carry and lift with ease



Compact – folds down to cabin luggage size



Portable – easy to take with you, anywhere



Robust – durable materials give peace of mind



Power
Chair

Compact, comfortable, quick and simple to fold... the Power Chair combines it all

- Easy to use, ultra-portable folding Power Chair
- Exceptionally lightweight at only 15kg*, yet remarkably robust, carrying passengers up to 135kg/21 Stone
- Folds with ease, fits into most car boots and is Airline-safe
- Ergonomically designed user controls, plus integrated attendant control for total flexibility
- On- or off-board charging options, and integral USB port for charging mobile phones
- The perfect travel companion!



The lighter way to empower your life

Enjoy newfound freedom and independence on your travel adventures

"With an eFOLDi you can just fold it up, and take it with you on buses, trains, or planes. Perfect for me to get out and about with my family and even go on holiday abroad!"

Anna, Reigate



*Excluding battery

Ultra-lightweight construction for strain-free lifting

Attendant Control comes as standard

Ergonomic, easy to use joystick controls

Ionised paintwork for exceptional durability

Electric and manual brake systems for maximum safety/control



Folds in seconds

High grip tyres for consistent traction

Anti-tippers for safety and reassurance

Easy to lift, store and transport



Rapid-charge, high performance lithium batteries



Specifications*

Maximum Speed	4mph / 6kph
Maximum Range	8 miles / 13km
Maximum User Weight	21 Stone / 135kg / 297lbs
Climbing Ability	1:10 Gradient
Folded Dimensions	22 x 13 x 36 inches / 58 x 34 x 92cm
Weight (excl. battery)	33lbs / 15kg
Weight (incl. battery)	37.5lbs/17kg
Type of Use	Pavement Only
Braking System	Electric and Manual Brake System
Airline friendly	Yes

* Specifications may vary due to manufacturing tolerances



Power Chair

...and enjoy being you.

"Being so portable and lightweight, it gave us so much more freedom to travel and explore."

Fraser Edward



Scan to see the Power Chair in action

0800 999 8123

sales@efoldi.com www.efoldi.com

Your local eFOLDi retailer:



SunTech UK Limited t/a eFOLDi, 25 Ormside Way, Holmethorpe Industrial Estate, Redhill, Surrey RH1 2LW.

eFOLDi